




**Universitätsklinikum
Tübingen**

MAX PLANCK INSTITUTE
FOR BIOLOGICAL CYBERNETICS



FASTPTX HANDS-ON SESSION

Dario Bosch^{1,2,3}

1 Core Facility MRT, Medical Faculty, University of Tübingen, Germany

**2 Department for Biomedical Magnetic Resonance, University of Tübingen,
Tübingen, Germany**

**3 High-Field MR Center, Max Planck Institute for Biological Cybernetics,
Tübingen, Germany**


ESMRMB Annual Meeting 2025, Marseille




PARALLEL TRANSMISSION IN ULTRAHIGH FIELD MRI

MRI at ultrahigh fields:

- Inhomogeneous B_1^+
- Highly localized tissue heating \rightarrow SAR
- Parallel Transmit Systems: many RF (B_1^+) transmitters



CP-mode excitation FA



pTx excitation FA

Very expensive high-end systems


- 9.4T (400 MHz) whole-body magnet
- 16 Tx channels
- Whole-body gradient coil
- In-house RF coil development



THE FastPtx OPTIMIZATION TOOLBOX

- For each time point: 16 complex RF channels + 3 gradient channels = 35 parameters

- Python → PyTorch
- automatic differentiation
→ efficient
- Supports STA and Bloch simulations
- Can run on CPU or GPU
- Supports individual or universal pulses



FastPtx

Freedom to change all parameters



Use the Pulse in Any Sequence



C²P Exchange



Views



Overview

Details

Create Announcement

Select Providing Type:
IDEA C2P

3251

Pulseseq on Philips MRI Systems: *Unlocking and Validating Open-Source Sequences*

Thomas Roos¹, Edwin Versteeg¹, Mark Gosselink¹, Hans Hoogduin¹, Dennis Klomp¹, Jeroen Siero¹, and Jannie Wijnen¹
¹Department of Radiology, High Field MRI group, University Medical Center Utrecht, Utrecht, Netherlands

Synopsis

Keywords: Pulse Sequence Design, Pulse Sequence Design, Pulseseq, reproducible, open, open-source, UHF, 7T, Philips

Motivation: Conventional MRI vendor-specific sequence development limits research transferability and reproducibility. Harmonized frameworks, like Pulseseq, overcome these limitations - but Pulseseq was not supported yet on Philips MRI scanners.

Goal(s): To develop a fully-compliant Pulseseq interpreter for Philips MRI systems, facilitating the use of open-source sequences previously incompatible with these platforms.

Approach: Custom adaptation of Pulseseq for the Philips R5.4 platform, followed by validation with standard imaging sequences using a field camera and phantom/in vivo scans, and comparison with native sequences.

Results: Demonstration of the first Philips-compatible Pulseseq interpreter, evidenced by successful scans at 7T, marking a leap in cross-vendor research capabilities.

Impact: The successful adaptation of Pulseseq sequences to Philips MRI system enables researchers to deploy and disseminate advanced MRI techniques universally, fostering cross-vendor collaboration and accelerating the evolution of MRI technology.



Liberating pTx from Vendor Lock-in: Open-Source Cross-Vendor Parallel Transmit MRI Sequences by extending pTx-Pulseseq to Siemens UHF Scanners

Dario Bosch^{1,2,3}, Thomas Roos⁴, and Klaus Scheffler^{1,2}

¹ High-Field MR Center, Max Planck Institute for Biological Cybernetics, Tübingen, Germany ² Department for Biomedical Magnetic Resonance, University of Tübingen, Tübingen, Germany ³ Core Facility for Magnetic Resonance Imaging, Medical Faculty, University of Tübingen, Germany ⁴ Department of Radiology, High Field MRI group, University Medical Center Utrecht, Utrecht, Netherlands

Introduction

MRI Pulse Sequence source code is **vendor-specific** and **sensitive to software updates** within one generation of MRI scanners. Pulseseq[1] is an approach for **vendor-agnostic open-source sequence programming**. Parallel Transmit (pTx) RF-Pulses are an important technique in ultra-high field MRI and were **not supported** by Pulseseq in the past, so **no open-source pTx**

Methods

Implementation

The **Pulseseq v1.4.3** interpreter for **Siemens MRI** scanners was modified to handle **pTx pulses**. The waveforms of all channels are concatenated and **stored as one long waveform** (Roos et al.[2]). Modern Siemens MRI systems require the **pTx pulses in .ini** files. The interpreter **converts pulses to .ini**

Pulseseq interpreter with pTx support

UHF

VE12U-AP04 F9, VE12U-SP01, VE12U-SP01_pTx, XA60A,

Dario Bosch

Edit




Fig. 1: The .seq file is used to create a .seq file, which converts the pTx pulse out of the pTx pulse.




Fig. 2: An open-source pTx workflow, complete from B1 mapping to image acquisition. With pTx support in Pulseseq, all sequences can be open-source.





Fig. 3: Individual channel GRE images in a phantom, recorded with the Pulseseq interpreter (left) and a vendor-provided sequence (right)



Funding by ERC Advanced Grant "SpreadMRI" No. 834940 is gratefully acknowledged.


ISMRM Ultra-High Field MR & Brain Function Workshop, Annapolis, 31.03-02.04.2025

CP

S1

S2

S3




TP

CP

S1

S2

S3



TP

- Toolbox and example data publicly available
- Bosch & Scheffler Magn Reson Mater Phy 2024.
<https://doi.org/10.1007/s10334-023-01134-7>
- github.com/dabosch/FastPtx



https://colab.research.google.com/github/dabosch/FastPtx/blob/main/Colab_Demo.ipynb

